



Airox™ Legendair™ Ventilator A mixed portable ventilator, integrating pressure and volume modes for the treatment of most acute and chronic respiratory pathologies, including COPD and neuromuscular diseases.

#### LEGENDAIR APPLICATIONS

- Home care
- Adults and infants weighing ≥5 kg
- Pressure and volume ventilation
- Noninvasive and invasive ventilation

#### **OXYGEN ENRICHMENT AND MONITORING**

The Legendair ventilator is equipped to receive and monitor oxygen supplied from an external low pressure source. The  ${\rm FiO_2}$  concentration can reach up to 50% of oxygen, and minimum and maximum alarms can be set.

### **V**, TARGET FEATURE

 $V_{_{\mathrm{T}}}$  Target: The  $V_{_{\mathrm{T}}}$  Target feature allows you to set a volume target within pressure control mode, thanks to automatic pressure adjustment. As patients sleep and move about,  $V_{_{\mathrm{T}}}$  Target works to maintain the established tidal volume setting (when in pressure mode) to improve alveolar ventilation optimization, potentially reducing nocturnal hypoventilation.

#### **COMFORT AND CONVENIENCE**

- Responsive flow trigger allows the ventilator to quickly adapt to the patient's breathing
- Small turbine promotes optimal synchronization of the patient to the ventilator
- $V_{_{\mathrm{T}}}$  Target: Ability to set target volume with maximum pressure
- Precise flow and pressure sensors with adjustable I and E triggers in pressure modes
- Quiet device: less than 30 dBA at 10 mbar

#### **PORTABILITY**


The Legendair ventilator integrates lithium-ion batteries that provide up to 10 hours of portable operation. It is also equipped with a 24 V DC input (Open Pack™ connection) for additional battery life and mobility.

#### ALARMS AND MONITORING MENU

The Legendair ventilator offers comprehensive alarm screens:

- Settings of minimum and maximum alarm values
- Display of the real measured ventilation values in the active cycles
- Display of the active alarm in a specific window
- Reading of the last nine stored alarm messages with date and hour stamp
- All alarms are stored in the memory and can be retrieved later






# **SIMPLIFIED MAINTENANCE**

- Quick and easy access to all technical components
- Micro-turbine with no maintenance for 15,000 hours
- RS 232 port for software upgrading
- Maintenance menu allows several diagnostics and sensor calibration
- Event memory allows viewing of ventilation alarms, technical events and other occurrences on the display
- Display shows patient and device hour meters




# **SETTINGS IN PSV/PSV ST MODE**

I <sub>PAP</sub> (inspiratory pressure)	5 - 40
E <sub>PAP</sub> (expiratory pressure)	0 - 20
Inspiratory ramp (Rise time)	4 Levels
Inspiratory trigger	5 Levels
Exhalation trigger	- 15 to - 75% of decrease from peak flow
Minimum rate	4 – 40
Apnea Time	3 – 30 s
V <sub>,</sub> Target	50 – 1400 mL
Maximum pressure (to reach V <sub>T</sub> Target)	8 to 55 mbar (limited to $+$ 3 to $+$ 15 of $I_{PAP}$ )
Monitoring	Vti, Vte, Rate, I/T Ratio, ${\rm FiO_2}$ , ${\rm I_{\rm PAP}}$ ${\rm E_{\rm PAP}}$ and patient triggering




### **SETTINGS IN PCV/PACV MODE**

I <sub>PAP</sub> (inspiratory pressure)	5 - 40
E <sub>PAP</sub> (expiratory pressure)	0 - 20
Inspiratory ramp (Rise time)	4 Levels
Minimum rate	5 – 60
I/T Ratio	50 to 25%
Inspiratory trigger	5 Levels
V <sub>,</sub> Target	50 – 1400 mL
Maximum pressure (to reach V <sub>T</sub> Target)	8 to 55 mbar (limited to $+$ 3 to $+$ 15 of $I_{PAP}$ )
Monitoring	Vti, Vte, Rate, I/T Ratio, ${\rm FiO_2}$ , ${\rm I_{\rm PAP}}$ , ${\rm E_{\rm PAP}}$ and patient triggering




## **SETTINGS IN CV/ACV MODE**

Tidal Volume	50 – 1400 mL
E <sub>PAP</sub> (expiratory pressure)	0 - 20
Flow ramp	Rectangular, Decelerated, Sinusoidal
Minimum rate	5 – 60
I/T ratio	50 to 25%
Inspiratory trigger	5 Levels
Sigh	Yes / No $V_{\tau}$ : from 1,0 to 2,0 of set $V_{\tau}$ Rate : 50 to 250
Monitoring	I <sub>PAP</sub> ,E <sub>PAP</sub> Vti, Vte, I/T, FiO <sub>2</sub> and patient triggering



### **SETTINGS IN SIMV MODE**

Tidal Volume	50 – 1400 mL
I <sub>PAP</sub> (inspiratory pressure)	5 - 40
E <sub>PAP</sub> (expiratory pressure)	0 - 20
SIMV Rate	2 - 20
I/T Ratio	50 to 25%
Inspiratory trigger	5 Levels
Security rate	4 - 40
Apnea Time	3 – 30 s
Monitoring	I <sub>PAP</sub> , E <sub>PAP</sub> , Vti, Vte, I/T, FiO <sub>2</sub> and patient triggering



#### **USER-FRIENDLY DESIGN**

- Large LCD screen
- All parameters displayed on a single screen
- Dedicated menu for alarms and patient monitoring
- · Ventilation mode can be changed during ventilation
- Tutorial messages for the user
- Parameters can be locked

# EXHALATION $V_{\downarrow}$ MONITORING

The Legendair ventilator offers the option to integrate exhaled tidal volume monitoring.

- Compatible with a variety of exhalation valves
- Exhalation kit is configured as a "mono-block" for easy installation and disinfection
- Exhaled tidal volume is shown on the display

#### **ERGONOMICS**

The Legendair ventilator was thoughtfully designed with numerous features that make it easy to handle and use.

- Large integrated handle
- Small size
- Lightweight: only 4.5 kg, including the internal battery
- Clear, easy-to-read display
- Patient circuits are quick and easy to install






### **TECHNICAL SPECIFICATIONS CHART**

Medical device class	II B
Weight	4.5 Kg
Dimensions	H 154 x L 235 x W 315 mm
Power supply	115/230 V ± 10% – 50/60 Hz
Isolation class	II .
Consumption	80 VA nominal and 90 VA Maxi
External battery	24 V DC
Internal battery	Lithium Ion quick charge.
6 to 10 hours autonomy	

# **ACCESSORIES AND PRODUCT REFERENCES CHART**

REFERENCES	DESCRIPTION
4095700	<b>Legendair Single Branch</b> (Equipped with: power supply cable, user's manual, carrying bag, single use patient circuit with exhalation valve)
4095701 (UK)	<b>Legendair Single Branch</b> (Equipped with: power supply cable, user's manual, carrying bag, single use patient circuit with exhalation valve)
4095900	<b>Legendair Double Branch</b> (Equipped with: power supply cable, user's manual, carrying bag, single use patient circuit with exhalation valve)
4095901 (UK)	<b>Legendair Double Branch</b> (Equipped with: power supply cable, user's manual, carrying bag, double branch patient circuit with exhalation valve)
	ACCESSORIES
4096099	External remote alarm (5 m)
3814199	FiO <sub>2</sub> kit (oxygen cell, T adaptator, connection cable, deflector)
3823099	Exhalation block
3809099	Carrying bag
2962799	Oxygen inlet connector (pack of 10)
3810899	24VDC external cable
2961999	Cable 9MF (RS232)
2967299	Dual bag
	CIRCUITS
5093600	Adult single limb circuit with exhalation valve 180 cm – single use
5094000	Adult double limb circuit with exhalation valve 180 cm – single use
5093500	Pediatric single limb circuit with exhalation valve 180 cm – single use
5093900	Pediatric double limb circuit with exhalation valve 180 cm – single use
	DISPOSABLES
2964299	Oxygen cell (for FiO <sub>2</sub> measurement kit)
2963399	Air inlet combi filter (pack of 6)





COVIDIEN, COVIDIEN with logo and "positive results for life" are trademarks of Covidien AG. TM marked brands are owned by Airox SA. © 2008 Covidien. - C-VE-Legendair/GB - 06/2009